

Pres-N-Lok Demountable Cell



ICL's Pres-N-Lok Cell enables the user to analyze mulls and smears without the use of clumsy springs, nuts, and bolts. The cell is assembled by placing the sample on one (1) of two (2) 25mm x 4mm crystal windows separated by a PTFE spacer, dropping them into the bore of the inner barrel, and pressing in the outer barrel. The cell consists of a back plate of anodized aluminum, a self locking inner barrel with a Viton® O-ring over which the outer barrel slips on and off with ease, and assorted PTFE spacers.

0006-2467, Pres-N-Lok Cell Complete with Assorted Spacers, but without Windows

Designed for optimum performance, ICL's Pres-N-Lok cell enables the user to analyze mulls and smears without the use of clumsy springs, nuts, bolts etc. Assembly of the cell requires only two crystal windows 25mm diameter by 4mm thick, one PTFE spacer and the user's sample. The cell consists of a back plate of anodized aluminum, a self locking barrel with a Viton® O-ring which slips the cover on and off with ease and assorted PTFE spacers. The cell pathlength can be varied through the use of different spacer thicknesses.

Additional Information

ICL's Pres-N-Lok Cell enables the user to analyze mulls and smears without the use of clumsy springs, nuts, and bolts. The cell is assembled by placing the sample on one (1) of two (2) 25mm x 4mm crystal windows separated by a PTFE spacer, dropping them into the bore of the inner barrel, and pressing in the outer barrel.

Windows : (need to order 2 windows , for ZnSe please contact us)

0002A-67	NaCl Polished Disc- 25mmx4mm	25 x 4 mm	NaCl
0002B-136	KCl Polished Disc- 25mmx4mm	25 x 4 mm	KCl
0002C-74	KBr Polished Disc- 25mmx4mm	25 x 4 mm	KBr
0002D-181	CaF2 Polished Disc- 25mmx4mm	25 x 4 mm	CaF2
0002E-176	BaF2 Polished Disc- 25mmx4mm	25 x 4 mm	BaF2
0002F-137	CsI Polished Disc- 25mmx4mm	25 x 4 mm	CSI
0002G-73	KRS-5 Polished Disc- 25mmx4mm	25 x 4 mm	KRS-5
0002H-321	Infrasil Quartz Polished Disc- 25mmx4mm	25 x 4 mm	Infrasil Quartz
0002K-1237	AgCl Polished Disc- 25mmx4mm	25 x 4 mm	AgCl

0006-4407, Pres-N-Lok Cell O-Rings- 25mm, pk/12



JASCO France
17 Rue des Cerisiers
91090 Lisses - France

Tél. 01 64 97 09 60 Fax: 01 85 09 90 76
email: jascofrance@jascofrance.fr