

Model G-2SS 10cm Stainless Steel Gas Cell



The G-2SS stainless steel gas cell is useful for qualitative analysis of pure gases and mixtures where one or more of the constituent gases are present in concentrations of at least 5 percent. It can be used in situations that make use of a glass cell impractical, such as applications involving high pressure or gases that etch a glass cell body.

The only functional constraint on the pressure that can be applied to this cell is the window material selected. The clear aperture of the cell is 39mm and it uses 47mm x 6mm windows. Like the G-2, this cell uses a universal 2" x 3" slide mount that fits all spectrophotometers. However, it is also available base plate mounted (See Fig. 1). There are inlet and outlet ports equipped with vacuum valves and hose barbs that are configured to also allow flowing gas through the cell. An 1/4 NPT port in the side of the cell allows the addition of a pressure gauge, a thermocouple and other probes. The cell uses nitrile "O" rings and PTFE gaskets. KBr and NaCl are the basic window material, but other materials are available upon request. Can be heated.



JASCO France
17 Rue des Cerisiers
91090 Lisses - France

Tél. 01 64 97 09 60 Fax: 01 85 09 90 76
email: jascofrance@jascofrance.fr

0008-5160 Model G-2SS 10cm Stainless Steel Gas Cell- NaCl

0008-5161 Model G-2SS 10cm Stainless Steel Gas Cell- KBr

0008-5162 Cylindrical Cell Holder for Model G-2SS 10cm Stainless Steel Gas Cell, Model G-3 PTFE 10cm PTFE Gas Cell and 10cm Vapor Gas Multi-Sampler

0008-5584 Optional Pressure Gauge for Model G-2SS 10cm Stainless Steel Gas Cell

0008-5804 PTFE Gaskets for Model G-2SS 10cm Stainless Steel Gas Cell- Pk/4

0008-5321 Nitrile 'O' Rings for Model G-2SS 10cm Stainless Steel Gas Cell- Pk/4

0008-5805 PTFE 'O' Rings for Model G-2SS 10cm Stainless Steel Gas Cell- Pk/4

0008-7329 Heater Jacket for 10cm PTFE Gas Cell and Model G-2SS 10cm Stainless Steel Gas Cell- 115v

0008-9404 Base Plate Mounting Option for G-2SS Stainless Steel Gas Cell



JASCO France
17 Rue des Cerisiers
91090 Lisses - France

Tél. 01 64 97 09 60 Fax: 01 85 09 90 76
email: jascofrance@jascofrance.fr

