





CATALOGUE VIALS 2021-2022

SOLUTIONS POUR LA CHROMATOGRAPHIE



JASCO

UltraVial

Resealable box helps to prevent contamination during use.
 Constructed with PTFE/ Silicone septa for contamination free sampling.
 Ultra Clean PTFE/ Silicone septa is optimized for gas chromatography and mass spectrometry.

> UltraVial - 2mL Vials & Closures



Ref.No.	JAS-01170	JAS-01172
Description	2mL Clear Glass 12x32mm Flat Base 9-425 Screw Thread Vial with Label. Blue 9-425 Open Top Ribbed Screw Cap with 9mm Red PTFE/White Silicone Septa 1mm Thick (UltraClean).	2mL Amber Glass 12x32mm Flat Base 9-425 Screw Thread Vial with Label. Blue 9-425 Open Top Ribbed Screw Cap with 9mm Red PTFE/White Silicone Septa 1mm Thick (UltraClean).

Vial Cap Septa are pre-assembled together. Blue Kit Packing. 100sets/pk



Ref.No.	JAS-01017	JAS-01018	JAS-01019
Description	Black 8-425 Open Top Screw Cap with 8mm Red PTFE/White Silicone Septa 1.5mm Thick (UltraClean).	Blue 11mm Open Top Snap Cap with Red PTFE/White Silicone Septa 1mm Thick (UltraClean).	Blue 9-425 Open Top Ribbed Screw Cap with 9mm Red PTFE/White Silicone Septa 1mm Thick (UltraClean).

100pcs/pk



Ref.No.	JAS-01020	JAS-01021	JAS-01028
Description	Blue 9-425 Open Top Ribbed Screw Cap with 9mm Red PTFE/White Silicone Septa 1mm Thick PRE-SLIT (UltraClean).	Red 9-425 Open Top Smooth Screw Cap with 9mm Red PTFE/White Silicone 1mm Thick (UltraClean).	Blue 9-425 Open Top Ribbed Screw Cap with 9mm White PTFE/Red Silicone Septa 1mm Thick (UltraClean).

100pcs/pk



Ref.No.	JAS-01037	JAS-01056
Description	Black 8-425 Open Top Screw Cap with 8mm Blue PTFE/White Silicone Septa 1.5mm Thick (UltraClean).	Blue 9-425 Open Top Ribbed Screw Cap with 9mm Blue PTFE/White Silicone Septa 1mm Thick PRE-SLIT (UltraClean).

100pcs/pk



> UltraVial - ND18



Ref.No.	JAS-00974	JAS-00975
Description	18mm Silver Color Open Top Metal Cap (8mm hole) with 17.5mm White PTFE/Blue Transparent Silicone Septa 1.5mm Thick (UltraClean).	18mm Silver Color Open Top Metal Cap (8mm hole) with 17.5mm Natural PTFE/Blue Transparent Silicone Septa 1.5mm Thick (UltraClean).
	100pcs/pk.	

> UltraVial - ND20



Ref.No.	JAS-00927	JAS-00972	JAS-00973
Description	20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm White PTFE/Blue Silicone Septa 3mm Thick (UltraClean).	20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm White PTFE/Blue Transparent Silicone Septa 3mm Thick (UltraClean).	20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm Natural PTFE/Blue Transparent Silicone Septa 3mm Thick (UltraClean).
	100pcs/pk		

> UltraVial - ND24



Ref.No.	JAS-00316	JAS-01082	JAS-01189
Description	22mm Blue PTFE/White Silicone 3.0mm thick Septa (Ultra Clean).	22mm Natural PTFE/White Silicone 3.0mm thick Septa (Ultra Clean).	24-400 White Open Top PP Screw Cap with 22mm Natural PTFE/ White Silicone 3.0mm Thick Septa (EPA Quality).
	100pcs/pk		



EPA/TOC Vials & Closures

If you are doing any sampling, testing, or validation, you need to make sure the results were not skewed by a contaminant that may have come from the container.

We offers a range of certified EPA Screw Vial Convenience Kits which are suitable on all TOC manufacturer's instruments.

All Certified to contain <10ppb Carbon, giving you confidence in the quality of your TOC measurements. Teflon lined caps and protective dust covers.

Superior vial-to-vial consistency.

Capped with teflon lined silicone septa and dust cover.

Widely apply to air, water, soil analysis test, also pesticide residue, pharmaceutical manufacturing and forensic test.

> UltraVial - EPA/TOC



Ref.No.	JAS-00413	JAS-00331	JAS-00338/JAS-00339
Description	20mL 27.5x57mm Clear Glass EPA /TOC Vial 24-400 White Open Top PP Screw Cap with 22mm Natural PTFE/White Silicone 3.0mm thick Septa (EPA Quality).	20mL 27.5x57mm Amber Glass EPA/TOC Vial 24-400 White Open Top PP Screw Cap with 22mm Natural PTFE/White Silicone 3.0mm thick Septa (EPA Quality).	20mL 27.5x57mm Clear/Amber Glass EPA/TOC Vial 24-400 White Open Top PP Screw Cap with 22mm Natural PTFE/White Silicone 3.0mm thick Septa (EPA Quality) and Dust Cover.
pre-assembled together. 100pcs/pk.			



Ref.No.	JAS-00334	JAS-00335	JAS-00342/JAS-00343
Description	40mL 27.5x95mm Clear Glass EPA/TOC Vial 24-400 White Open Top PP Screw Cap with 22mm Natural PTFE/White Silicone 3.0mm thick Septa (EPA Quality).	40mL 27.5x95mm Amber Glass EPA/TOC Vial 24-400 White Open Top PP Screw Cap with 22mm Natural PTFE/White Silicone 3.0mm thick Septa (EPA Quality).	40mL 27.5x95mm Clear/Amber Glass EPA/TOC Vial 24-400 White Open Top PP Screw Cap with 22mm Natural PTFE/White Silicone 3.0mm thick Septa (EPA Quality) and Dust Cover.
pre-assembled together. 72pcs/pk.			



Ref.No.	JAS-00336	JAS-00337	Blue Kit Packing
Description	60mL 27.5x140mm Clear Glass EPA/TOC Vial 24-400 White Open Top PP Screw Cap with 22mm Natural PTFE/White Silicone 3.0mm thick Septa (EPA Quality).	60mL 27.5x140mm Amber Glass EPA VOA VOC TOC Vial 24-400 White Open Top PP Screw Cap with 22mm Natural PTFE/White Silicone 3.0mm thick Septa (EPA Quality).	
pre-assembled together. 100pcs/pk.			

Bonded

Our bonding machines are fully automated to bond liners to caps using a mechanical bond without any adhesives or chemicals.

No glue contamination.

No additional assembly.

Prevent liner fallout.

Use with different polypropylene caps.

Available in 1mm-3mm standard thicknesses; customization is available.

Ideal for medical devices, diagnostics packaging and specialty chemicals.

>ND 9 Closures



Ref.No.	JAS-00144	C0001359
Description	Blue 9-425 Open Top Ribbed Screw Cap Bonded with 9mm Red PTFE/White Silicone Septa 1mm Thick.	Blue 9-425 Open Top Ribbed Screw Cap Bonded with 9mm Red PTFE/White Silicone Septa 1mm Thick PRE-SLIT.
		100pcs/pk.

>ND 24 Closures



Ref.No.	JAS-01414
Description	24-400 White Open Top PP Screw Cap Bonded with 22mm Natural PTFE/White Silicone Septa 3mm Thick.
	100pcs/pk.



AUTOSAMPLER VIALS & CLOSURES

Screw Neck ND8

Standard opening vials for GC and HPLC, used on many instruments.

To meet your micro sampling needs, we have a variety of micro volume inserts with capacities and designs that can transform our vials from full-capacity to limited volume with one simple step.

Standard opening requires Micro-insert with a diameter of 5mm.

Micro inserts with flat bottom, conical base, with plastic spring are available.

The conical design permits complete sample removal.



Ref.No.	JAS-00001	JAS-00002	JAS-00003	JAS-00004
Description	2mL Clear Glass 12x32mm Flat Base 8-425 Screw Thread Vial.	2mL Clear Glass 12x32mm Flat Base 8-425 Screw Thread Vial with Label.	2mL Amber Glass 12x32mm Flat Base 8-425 Screw Thread Vial.	2mL Amber Glass 12x32mm Flat Base 8-425 Screw Thread Vial with Label.
	100pcs/pk			

> 8mm Septa




Ref.No.	JAS-00436	JAS-00943	JAS-00321	JAS-00964	JAS-00445	JAS-00999
Description	8-425 Red PTFE/White Silicone Septa 1mm Thick.	8-425 Red PTFE/White Silicone Septa 1.3mm Thick.	8-425 Red PTFE/White Silicone Septa 1.5mm Thick.	8-425 Red PTFE/White Silicone Septa 1.5mm Thick Pre-slit.	8-425 White PTFE/Red Silicone Septa 1.5mm thick.	8-425 Blue PTFE/White Silicone Septa 1.5mm thick.
	100pcs/pk					

> ND8 Closures



Ref.No.	JAS-00931	JAS-00957	JAS-00133	JAS-00409
Description	Black 8-425 Open Top Screw Cap.	White 8-425 Open Top Screw Cap.	Black 8-425 Open Top Screw Cap with 8mm Red PTFE/White Silicone Septa 1.5mm Thick.	White 8-425 Open Top PP Cap with 8mm White PTFE/Red Silicone Septa 1.5mm Thick.
	100pcs/pk			

> 5mm Micro Insert



Ref.No.	JAS-00074	JAS-00075	JAS-00076
Description	5x29mm Insert Clear Glass Conical Base with Polyspring.	5x31mm Insert Clear Glass Conical Base.	5x31mm Insert Clear Glass Flat Base.
	100pcs/pk		

Short Neck ND9 (recommended for JASCO)

Wide opening 9-425 vials can be used on all common autosamplers, **perfectly compatible with JASCO** and many others

Vials with intergrated micro-insert are available in clear and amber glass.

Wide opening 12x32mm, 2ml Polypropylene vials have chemically resistance, which are used for atomic absorption, water and protein analysis, capillary electrophoresis, ion chromatography applications etc.

PTFE/silicone slit shape choices are Single, Y-shape, Cross available.

The polymer spring acts as a shock absorber that protects against breakage if the needle bottoms out.

Vials and closures are available for packing into a kit to reduce the risk of contamination of vials in laboratories.

A variety of septa colors and pre-assembled closures are available to offer as request.



Ref.No.	JAS-00006	JAS-00008	JAS-00010	JAS-00011
Description	2mL Clear Glass 12x32mm Flat Base 9-425 Screw Thread Vial.	2mL Clear Glass 12x32mm Flat Base 9-425 Screw Thread Vial with Label.	2mL Amber Glass 12x32mm Flat Base 9-425 Screw Thread Vial.	2mL Amber Glass 12x32mm Flat Base 9-425 Screw Thread Vial with Label.
100pcs/pk				



Ref.No.	JAS-00084	JAS-00410	JAS-00127
Description	0.3mL 9-425 PP Screw Thread Micro-Vial Transparent. 12x32mm.	2mL 9-425 Wide Opening PP Screw Vial with Graduations Transparent .	2mL 9-425 12x32mm Clear Screw Thread Vial Bottom with Integrated 0.2ml Glass Micro-insert.
100pcs/pk			



Ref.No.	JAS-00961	JAS-00962
Description	2mL Amber Glass 12x32mm Flat Base 9-425 Screw Thread Vial with Label. Blue 9-425 Open Top Ribbed Screw Cap with 9mm White PTFE/Red Silicone Septa 1mm Thick.	2mL Clear Glass 12x32mm Flat Base 9-425 Screw Thread Vial with Label. Blue 9-425 Open Top Ribbed Screw Cap with 9mm White PTFE/Red Silicone Septa 1mm Thick.
Kit Packing.100pcs/pk		

> 9mm Septa



Ref.No.	JAS-00402	JAS-00329	JAS-00439	JAS-00324	JAS-00431
Description	9-425 Red PTFE/White Silicone Septa 1mm Thick.	9-425 White PTFE/Red Silicone Septa 1mm Thick. .	9-425 Blue PTFE/White Silicone Septa 1mm Thick.	9-425 Red PTFE/White Silicone Septa 1mm Thick Pre-slit.	9-425 Blue PTFE/White Silicone Septa 1mm Thick Pre-slit.
	100pcs/pk				



Ref.No.	JAS-01133	JAS-00801	JAS-00460
Description	9-425 Red PTFE/White Silicone Septa 1.0mm Thick Pre-slit(Mercedes Shape).	9-425 Red PTFE/White Silicone Septa 1mm Thick Pre-slit(Cross Shape).	9-425 Natural PTFE/White Silicone Septa 1mm thick.
	100pcs/pk		

> ND9 Cap





Ref.No.	JAS-00421	JAS-00327	JAS-01132	JAS-01051	JAS-00856
Description	Blue 9-425 Open Top Smooth Screw Cap.	Blue 9-425 Open Top Ribbed Screw Cap.	Black 9-425 Open Top Smooth Screw Cap.	Black 9-425 Open Top Ribbed Screw Cap.	Gold Magnetic 9-425 Open Top Screw Cap.
	100pcs/pk				

> ND9 Closures



Ref.No.	JAS-00143	JAS-00146	JAS-00148	JAS-00149	JAS-00152
Description	Blue 9-425 Open Top Ribbed Screw Cap With 9mm Red PTFE/White Silicone Septa 1mm Thick	Blue 9-425 Open Top Ribbed Screw Cap With 9mm White PTFE/Red Silicone Septa 1mm Thick	Blue 9-425 Open Top Ribbed Screw Cap With 9mm Red PTFE/White Silicone Septa 1mm Thick Pre-slit	Blue 9-425 Open Top Ribbed Screw Cap With 9mm Blue PTFE/White Silicone Septa 1mm Thick Pre-slit	Blue 9-425 Open Top Ribbed Screw Cap With 9mm White PTFE/Red Silicone Septa 1mm Thick Pre-slit
	100pcs/pk				

> 6mm Micro Insert



Ref.No.	JAS-00070	JAS-00072	JAS-00073	JAS-00077
Description	6x29mm Insert Clear Glass Conical Base with Polyspring.	6x31mm Insert Clear Glass Conical Base.	6x31mm Insert Clear Glass Flat Base.	6x29mm PP Insert Conical Base.
	100pcs/pk			

Crimp Neck ND11

The Crimp Vials are used on instruments of JASCO and many others.

The wide opening crimp cap provides a larger target area for improved autosampler needle accuracy.

Vials with intergrated Micro-insert are also available in clear and amber glass.

PTFE/silicone, Natural Rubber/TEF septa, single slit septa to choose.

The 3-layer septa of Natural Rubber/Red Butyl/TEF combines the good physical properties of Nat. Rubber (Resealability) with good chemical properties of butyl. (See Customized service page)



Ref.No.	JAS-00014	JAS-00015	JAS-00016	JAS-00018	JAS-00129
Description	2mL Clear Glass Flat Base 11mm Crimp Vial Wide Opening.	2mL Clear Glass Flat Base 11mm Crimp Vial Wide Opening with Label.	2mL Amber Glass Flat Base 11mm Crimp Vial Wide Opening.	2mL Amber Glass Flat Base 11mm Crimp Vial Wide Opening with Label.	2mL 11mm Crimp Top Clear Vial Bottom with Integrated 0.2ml Glass Micro-insert.
	100pcs/pk				

ND 11 Closures



Ref.No.	JAS-00164	JAS-00165	JAS-00166	JAS-00167
Description	SILVER 11mm Open Top Crimp Cap with CLEAR TEF/ Red Silicone Septa 1mm Thick.	SILVER 11mm Open Top Crimp Cap with Red PTFE/ White Silicone Septa 1mm Thick.	SILVER 11mm Open Top Crimp Cap with Red PTFE/ White Silicone/Red PTFE Septa 1mm Thick.	SILVER 11mm Open Top Crimp Cap with Red Silicone Septa 1mm Thick.
	100pcs/pk			



Ref.No.	JAS-00168	JAS-00169★ Patent No.:CN107163311B	JAS-00980
Description	SILVER 11mm Open Top Crimp Cap with White PTFE/ Red Silicone Septa 1mm Thick.	SILVER 11mm Open Top Crimp Cap with Transparent TEF/Natural Rubber red-orange Septa 1mm Thick.	SILVER 11mm Open Top Crimp Cap with Red PTFE/ White Silicone Septa 1mm Thick Pre-slit.
	100pcs/pk		

Snap Neck ND11

The vials are used on instruments of JASCO and many others.

Snap Ring finish eliminates the need for crimping or decapping.

The intergrated Micro-insert are available for in clear and amber glass, micro-inserts can be delivered pre-assembled in vials.

Snap ring caps are also available in a soft quality with transparent or blue caps, soft caps are more convenient in handling (pushing on/removal).

The septa shapes have Single and Snow type to choose.



Ref.No.	JAS-00020	JAS-00021	JAS-00022	JAS-00023	JAS-00130
Description	2mL Clear Glass Flat Base 11mm Snap Vial Wide Opening.	2mL Clear Glass Flat Base 11mm Snap Vial Wide Opening with Label.	2mL Amber Glass Flat Base 11mm Snap Vial Wide Opening.	2mL Amber Glass Flat Base 11mm Snap Vial Wide Opening with Label.	2mL 11mm Snap Top Clear Vial Bottom with Integrated 0.2ml Glass Micro-insert.
	100pcs/pk				

>ND 11 Closures



Ref.No.	JAS-00170	JAS-00171	JAS-00172	JAS-00173
Description	CLEAR 11mm Open Top Snap Cap with Red PTFE/White Silicone Septa 1mm Thick.	BLUE 11mm Open Top Snap Cap with Red PTFE/White Silicone Septa 1mm Thick.	RED 11mm Open Top Snap Cap with Red PTFE/White Silicone Septa 1mm Thick.	CLEAR 11mm Open Top Snap Cap with Red PTFE/White Silicone/Red PTFE Septa 1mm Thick.
	100pcs/pk			



Ref.No.	JAS-00174	JAS-00175	JAS-00176	JAS-00177
Description	CLEAR 11mm Open Top Snap Cap with Red PTFE/White Silicone Septa 1mm Thick Pre-slit.	CLEAR 11mm Open Top Snap Cap with Red PTFE/White Silicone Septa 1mm Thick Pre-slit (snow shape).	CLEAR 11mm Open Top Snap Cap with Transparent TEF/Natural Rubber red-orange Septa 1mm Thick.	CLEAR 11mm Open Top Snap Cap with White PTFE/Red Silicone Septa 1mm Thick.
	100pcs/pk			

Short Thread ND13

The vials are used on instruments of JASCO and many others.

4ML vials are widely used in compound storage as well as for Chromatography vials, available as closed or open top screw seals with 13-425 thread.



Ref.No.	JAS-00024	JAS-00025	JAS-00026	JAS-00027
Description	4mL Clear Glass 15x45mm Flat Base 13-425 Screw Thread Vial.	4mL Clear Glass 15x45mm Flat Base 13-425 Screw Thread Vial with Label.	4mL Amber Glass 15x45mm Flat Base 13-425 Screw Thread Vial.	4mL Amber Glass 15x45mm Flat Base 13-425 Screw Thread Vial with Label.

100pcs/pk

>ND 13 Septa



Ref.No.	JAS-00322	JAS-00408	JAS-00863	JAS-00935	JAS-01134
Description	13mm Red PTFE/White Silicone Septa 1mm Thick.	13mm Red PTFE/White Silicone Septa 1mm Thick Pre-slit.	13mm White PTFE/Red Silicone Septa 1mm Thick.	13mm Natural PTFE/White Silicone Septa 1mm thick.	13mm Pure PTFE Septa 1.0mm Thick.

100pcs/pk

>ND 13 Closures



Ref.No.	JAS-00179	JAS-00180	JAS-00181	JAS-00182
Description	Black PP 13-425 Open Top Screw Cap with Red PTFE/White Silicone Septa 1mm Thick.	Black PP 13-425 Open Top Screw Cap with Red PTFE/White Silicone Septa 1mm Thick Pre-slit.	Black PP 13-425 Open Top Screw Cap with White PTFE/Red Silicone Septa 1mm Thick.	Black PP 13-425 Open Top Screw Cap with Natural PTFE/White Silicone Septa 1mm Thick.

100pcs/pk

Shell Vial

The vial with soft plug is used on JASCO and many others, the star-shaped plug is for easy penetration.

The normal size is 8x40mm, and 8x30mm is also available.



Ref.No.	JAS-00078	JAS-00080
Description	1mL 8x40mm Flat Base Clear Glass Shell Vial with Plug.	1mL 8x30mm Flat Base Clear Glass Shell Vial with Plug.
200pcs/pk.		

High Recovery Vial

High recovery vials combine large sample capacities with near full sample extraction.

Excellent for Compound Storage.


Smoothly sloped inner conical surface ensures a high level of solid/liquid recovery.



Ref.No.	JAS-01211	JAS-01212
Description	1.5mL High Recovery Vial Clear Glass 12x32mm Flat Base 9-425 Screw Thread.	1.5mL High Recovery Vial Amber Glass 12x32mm Flat Base 9-425 Screw Thread.
100pcs/pk.		

Test Tube

Screw thread test tubes include PP caps with PTFE/silicone septa, which is made from Borosilicate glass tubing.



Ref.No.	JAS-01093	JAS-01138	JAS-00131
Description	15-415 Black Closed Top PP Cap with Blue PTFE/White Silicone Septa 1.5mm Thick.	15-415 White Closed Top PP Cap for Test Tube with 15mm Natural PTFE/White Silicone Septa 1.5mm Thick.	Test Tube Clear Glass 16x100mm Rounded Bottom Capacity to Neck 11mL. Chinesing Tubing.
250pcs/pk.			

Screw Neck ND10

Wide-opening vials provide an increased target area for easier sample preparation and to reduce the chance of bent or broken needles during sampling.

Designed to work in Agilent and other autosamplers.



Ref.No.	JAS-01191	JAS-01192	JAS-01210
Description	2mL Clear Glass 12x32mm Flat Base 10-425 Screw Thread Vial with Label.	2mL Amber Glass 12x32mm Flat Base 10-425 Screw Thread Vial with Label.	Black 10-425 Open Top Screw Cap with Red PTFE/White Silicone Septa 1mm Thick.
100pcs/pk			

HEADSPACE VIALS & CLOSURES



Screw Headspace ND18

Pre-assembling vial+closure service is available upon request.
The magnetic screw caps can be used for headspace as well as for SPME.
Type of Cap: Magnetic Precisionj screwcap, Closed Precisionj screwcap.



Ref.No.	JAS-00028	JAS-00029	JAS-00030	JAS-00031	JAS-00920	JAS-01139
Description	10mL Clear Glass Screw Headspace Vial 22.5x46mm.	10mL Amber Glass Screw Headspace Vial 22.5x46mm.	20mL Clear Glass Screw Headspace Vial 22.5x75mm.	20mL Amber Glass Screw Headspace Vial 22.5x75mm.	18mm Silver Color Open Top Metal Cap.	18mm Silver Color Closed Top Metal Cap.
	100pcs/pk					



Ref.No.	JAS-00208	JAS-00474	JAS-00400	JAS-00438
Description	18mm Silver Color Open Top Metal Cap (8mm hole) with 17.5mm Blue PTFE/White Silicone Septa 1.5mm Thick.	18mm Silver Color Closed Top Metal Cap with Blue PTFE/White Silicone Septa 1.5mm Thick.	Blue PTFE/White Silicone Septa 1.5mm thick. (used for 18mm magnetic screw cap).	Red PTFE/White Silicone Septa 1.5mm Thick.
	100pcs/pk			

>ND 18 Pre-assembled



Ref.No.	JAS-00977	JAS-00976
Description	10mL Clear Glass Screw Headspace Vial 22.5x46mm 18mm Silver Color Open Top Metal Cap (8mm hole) with 17.5mm Blue PTFE/White Silicone Septa 1.5mm Thick Vial Cap Septa are pre-assembled together.	20mL Clear Glass Screw Headspace Vial 22.5x75mm 18mm Silver Color Open Top Metal Cap (8mm hole) with 17.5mm Blue PTFE/White Silicone Septa 1.5mm Thick Vial Cap Septa are pre-assembled together.
	100pcs/pk	

Crimp Headspace ND20

GC Headspace technology for the analysis of volatile organic solids and liquid samples after gasification
The headspace volume have 6ml, 10ml, 20ml clear/amber to choose, the bottom have flat, rounded and rounded-flat.

Type of Cap:

Aluminum Cap, Bimetal cap, Magnetic crimp cap.

Type of Liner:

Butyl, Butyl/PTFE, Pharma-fix liner, PTFE/silicone, Butyl stopper.



Ref.No.	JAS-00032	JAS-00033	JAS-00034	JAS-00035	JAS-00036
Description	6mL Clear Glass 22x38mm. 20mm Beveled Edge. Flat Bottom. Short Neck. Crimp Headspace Vial.	10mL Clear Glass 22.5x46mm. 20mm Beveled Edge. Flat Bottom. Crimp Headspace Vial.	10mL Clear Glass 22.5x46mm. 20mm Beveled Edge. Rounded-flat Bottom. Crimp Headspace Vial.	10mL Amber Glass 22.5x46mm. 20mm Beveled Edge. Flat Bottom. Crimp Headspace Vial.	20mL Clear Glass 22.5x75mm. 20mm Beveled Edge. Flat Bottom. Crimp Headspace Vial.
	100pcs/pk				



Ref.No.	JAS-00038	JAS-00039	JAS-00040
Description	20mL Clear Glass 22.5x75mm. 20mm Beveled Edge. Rounded-Flat Bottom. Crimp Headspace Vial.	20mL Clear Glass 22.5x75mm. 20mm Beveled Edge. Rounded-Flat Bottom. Crimp Headspace Vial. Short Neck.	20mL Amber Glass 22.5x75mm. 20mm Beveled Edge. Rounded-Flat Bottom. Crimp Headspace Vial. Short Neck.
	100pcs/pk		



> ND 20 Closures



Ref.No.	JAS-00196	JAS-00197	JAS-00198
Description	20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm Natural PTFE/White Silicone Septa 3mm Thick.	20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm Natural PTFE/Blue Silicone Septa 3mm Thick.	20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm White PTFE/White Silicone Septa 3mm Thick.
	100pcs/pk		



Ref.No.	JAS-00199	JAS-00200	JAS-00201
Description	20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm White PTFE/Blue Silicone Septa 3mm Thick.	20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm Red PTFE/White Silicone Septa 3mm Thick.	20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm Natural PTFE/Natural Silicone Septa 3mm Thick.
	100pcs/pk		



Ref.No.	JAS-00204	JAS-00205	JAS-00206
Description	20mm Open Top Gold Magnetic Crimp Cap (10mm hole) with 20mm Natural PTFE/Blue Silicone Septa 3mm Thick	20mm Open Top Gold Magnetic Crimp Cap (10mm hole) with 20mm Natural PTFE/White Silicone Septa 3mm Thick.	20mm Open Top Bi-Metallic Crimp Cap (8mm hole) Blue & Silver with 20mm Natural PTFE/White Silicone Septa 3mm Thick.
	100pcs/pk		



Ref.No.	JAS-00207	JAS-00835	JAS-00855
Description	20mm Open Top Bi-Metallic Crimp Cap (8mm hole) Blue & Silver with 20mm Natural PTFE/Natural Silicone Septa 3mm Thick.	20mm Open Top Gold Magnetic Crimp Cap (10mm hole) with 20mm White PTFE/Blue Silicone Septa 3mm Thick.	20.2mm Open Top Aluminum Crimp Cap (10mm hole) with Moulded Grey PTFE/Butyl Septa one ring 3mm Thick.
	100pcs/pk		



Ref.No.	JAS-00951	JAS-00952	JAS-01012
Description	20mm Open Top Bi-Metallic Crimp Cap (8mm hole) Red & Silver with 20mm Natural PTFE/Natural Silicone Septa 3mm Thick.	20mm Open Top Bi-Metallic Crimp Cap (8mm hole) Gold & Silver with 20mm Natural PTFE/Natural Silicone Septa 3mm Thick.	20.2mm Open Top Aluminum Crimp Cap (10mm hole) with Moulded Pharma-fix-septa Grey PTFE/Butyl 3mm thick.
	100pcs/pk		

> 20mm Butyl Septa



Ref.No.	JAS-00407	JAS-00433	JAS-00448	JAS-00449	JAS-00450	JAS-00423
Description	20mm Moulded Grey PTFE/Butyl Septa one ring 3mm thick.	20mm Moulded Grey Butyl Injection Stopper.	20mm Moulded Grey Butyl Septa one ring. 3mm thick.	20mm Moulded Grey Butyl Septa both rings.	20mm Grey PTFE/Butyl without ring. 3mm thick.	20mm Moulded Pharma-fix-septa Grey PTFE/Butyl 3mm thick.
	100pcs/pk					

Sample Vials & Closures

> Sample Vials & Closures

Screw vials for storage purpose, vials with different volumes are available.

Vial Color: Clear, Amber.

Cap Color: White, Black, Red, Blue, Yellow, Green, etc.

Septa material: PE, PTFE/Silicone, Silicone only, etc.

Pre-assemble service is available.



Ref.No.	JAS-00006	JAS-00041	JAS-00024	JAS-00044	JAS-00048	JAS-00049
Description	2mL Clear Glass Sample Vial 12x32mm 9-425 Screw Thread.	2mL Amber Glass Sample Vial 12x32mm 9-425 Screw Thread.	4mL Clear Glass Sample Vial 15x45mm 13-425 Screw Thread.	4mL Amber Glass Sample Vial 15x45mm 13-425 Screw Thread.	7mL Clear Glass Sample Vial 22x40mm 18-400 Screw Thread.	7mL Amber Glass Sample Vial 22x40mm 18-400 Screw Thread.
100pcs/pk						



Ref.No.	JAS-00051	JAS-00052	JAS-00053	JAS-00054	JAS-00055	JAS-00056
Description	8mL Clear Glass Sample Vial 17x60mm 15-425 Screw Thread.	8mL Amber Glass Sample Vial 17x60mm 15-425 Screw Thread.	10mL Clear Glass Sample Vial 22x52mm 18-400 Screw Thread.	10mL Amber Glass Sample Vial 22x52mm 18-400 Screw Thread.	16mL Clear Glass Sample Vial 22x72.5mm 18-400 Screw Thread.	16mL Amber Glass Sample Vial 22x72.5mm 18-400 Screw Thread.
100pcs/pk						



Ref.No.	JAS-00057	JAS-00058	JAS-00059	JAS-00060	JAS-01222	JAS-01223
Description	20mL Clear Glass Sample Vial 27.5x57mm 24-400 Screw Thread.	20mL Amber Glass Sample Vial 27.5x57mm 24-400 Screw Thread.	30mL Clear Glass Sample Vial 27.5x84mm 24-400 Screw Thread.	30mL Amber Glass Sample Vial 27.5x84mm 24-400 Screw Thread.	22mL Clear Glass Sample Vial 23x85mm. 20-400 Screw Thread.	22mL Amber Glass Sample Vial 23x85mm 20-400 Screw Thread.
100pcs/pk						




Ref.No.	JAS-00061	JAS-00062	JAS-00065	JAS-00066
Description	40mL Clear Glass Sample Vial 27.5x95mm 24-400 Screw Thread.	40mL Amber Glass Sample Vial 27.5x95mm 24-400 Screw Thread.	60mL Clear Glass Sample Vial 27.5x140mm 24-400 Screw Thread.	60mL Amber Glass Sample Vial 27.5x140mm 24-400 Screw Thread.
	100pcs/pk			

> Storage Cap & Septa



Ref.No.	JAS-00453	JAS-00455	JAS-00864	JAS-01135	JAS-00858	JAS-00482
Description	9-425 Blue Smooth Closed PP Top Cap.	13-425 Black Closed Top PP Cap.	13-425 White Closed Top PP Cap.	13-425 Red Closed Top PP Cap.	13-425 Green Closed Top PP Cap.	15-425 Black Closed Top PP Cap.
	100pcs/pk					



Ref.No.	JAS-00988	JAS-00480	JAS-00456	JAS-01136	JAS-00859	JAS-00860
Description	15-425 White Closed Top PP Cap.	18-400 Black Closed Top PP Cap.	18-400 White Closed Top PP Cap.	20-400 Black Closed Top PP Cap.	20-400 White Closed Top PP Cap.	22-400 Black Closed Top PP Cap.
	100pcs/pk					



Ref.No.	JAS-00475	JAS-00481	JAS-00457	JAS-00402	JAS-00322
Description	22-400 White Closed Top PP Cap.	24-400 Black Closed Top PP Cap.	24-400 White Closed Top PP Cap.	9-425 Red PTFE / White Silicone Septa 1.0mm Thick.	13-425 Red PTFE / White Silicone Septa 1.0mm Thick.
	100pcs/pk				



Ref.No.	JAS-00458	JAS-00440	JAS-00958	JAS-00862	JAS-00443
Description	15-425 Natural PTFE/White Silicone Septa 1.0mm Thick.	18-400 Natural PTFE/White Silicone Septa 1.5mm Thick.	20-400 Natural PTFE/White Silicone Septa 2.0mm Thick.	22-400 Natural PTFE/White Silicone Septa 2.0mm Thick.	24-400 Natural PTFE/White Silicone Septa 1.5mm Thick.
	100pcs/pk				

24-400 Closures With Sample Vial

If you are doing any sampling, testing, or validation, you need to make sure the results were not skewed by a contaminant that may have come from the container.

Teflon lined caps and protective dust covers.

Superior vial-to-vial consistency.



Ref.No.	JAS-01144	JAS-01417
Description	20mL Clear Glass Sample Vial. 27.5x58mm with 24-400 Black Closed Top PP Cap with 22mm Natural PTFE/White Silicone Septa 1.5mm thick.	20mL Amber Glass Sample Vial. 27.5x57mm. with 24-400 White Closed Top PP Cap with 22mm Natural PTFE/White Silicone Septa 1.5mm thick.

Vial+Cap+Septa are pre-assembled together.100pcs/pk.



Ref.No.	JAS-01096	JAS-01125
Description	40mL Clear Glass Vial 24-400 White Open Top PP Screw Cap and Natural PTFE/White Silicone 3mm Thick Septa.	40mL Amber Glass Vial 24-400 White Open Top PP Screw Cap and Natural PTFE/White Silicone 3mm Thick Septa.



Vial + Cap + Septa will be pre-assembled together. 72pcs/pk.



Ref.No.	JAS-01206	JAS-01207
Description	60mL Clear Glass Vial 24-400 White Open Top PP Screw Cap and Natural PTFE/White Silicone 3mm Thick Septa.	60mL Amber Glass Vial 24-400 White Open Top PP Screw Cap and Natural PTFE/White Silicone 3mm Thick Septa.

Vial + Cap + Septa will be pre-assembled together. 100pcs/pk.

Sample Vial Kit Package





Ref.No.	JAS-01316	JAS-01317
Description	4mL Clear Glass Sample Vial with 13-425 Red Closed Top PP Screw Cap and Red PTFE/White Silicone 1mm Thick Septa.	4mL Amber Glass Sample Vial with 13-425 Red Closed Top PP Screw Cap and Natural PTFE/White Silicone 1mm Thick Septa.
Vial + Cap + Septa are pre-assembled together. Packed with dividers. 144pcs/pk.		



Ref.No.	JAS-01327	JAS-01319
Description	8mL Clear Glass Sample Vial 17x60mm with 15-425 White Closed Top PP Cap with Natural PTFE/White Silicone Septa 1.0mm thick.	8mL Amber Glass Sample Vial 17x60mm with 15-425 White Closed Top PP Cap with Natural PTFE/White Silicone Septa 1.0mm thick.
Vial + Cap + Septa are pre-assembled together. Packed with dividers. 144pcs/pk.		




Ref.No.	JAS-01328	JAS-01321
Description	20mL Clear Glass Sample Vial. 27.5x58mm with 24-400 White Closed Top PP Cap with 22mm Natural PTFE/White Silicone Septa 1.5mm thick.	20mL Amber Glass Sample Vial. 27.5x58mm with 24-400 White Closed Top PP Cap with 22mm Natural PTFE/White Silicone Septa 1.5mm thick.
Vial + Cap + Septa are pre-assembled together. Packed with dividers. 144pcs/pk.		



Ref.No.	JAS-00802	JAS-00404
Description	20mL Clear Glass Vial. 27.5x58mm with 22-400 White Closed Top PP Cap and Silver Liner.	20mL Amber Glass Vial. 27.5x58mm with 22-400 White Closed Top PP Cap and Natural PTFE / White Silicone Septa.
Vial + Cap + Septa are pre-assembled together. Packed with dividers. 100pcs/pk		

Téléchargez nos brochures Chromatographie et Filtres Seringues



Télécharger



Télécharger

Contactez-nous pour en savoir plus sur nos produits ou consultez notre site Internet



JASCO FRANCE SAS
17 rue des Cerisiers
91090 LISSES

01 64 97 09 60
jascofrance@jascofrance.fr
www.jascofrance.fr